

Educational Resource



Health
Hunter New England
Local Health District

Wound Management: Inadine

Relevant for use at clinical sites where the product is used

Target: Clinicians **My Health Learning:** N/A

Purpose: Information/ demonstration of application and use of the product

Description: This resource provides information/ demonstration on the use of the product using a presentation. These can be used a part of a broader training program for wound management or as a standalone resource to inform the use of this specific product.

Learning outcomes: Identify indication/contraindications of Inadine
Describe how to apply Inadine

Keywords: Inadine, iodine dressing

Related Legislation, Australian Standard, NSW Ministry of Health Policy Directive or Guideline, National Safety and Quality Health Service Standard (NSQHSS) and/or other, HNE Health Document, Professional Guideline, Code of Practice or Ethics:

NSQHS Standards Criterion and/or other: 1, 3, 8

[NSW Health Policy Directive Clinical Procedure Safety PD 2017_032](#)

[NSW Ministry of Health Policy PD 2005_406 Consent to Medical Treatment](#)

[NSW Ministry of Health Policy Directive PD 2017_032 Infection Control Policy](#)

[Wounds Australia 2016 Standards for Wound Management Prevention and Management 3rd Ed](#)

Contact: Jayne Campbell jayne.campbell@hnehealth.nsw.gov.au

Issued: August 2015 **Review:** August 2018

How To Use Inadine


[Double click to view the presentation](#)

INADINE®

What is it?
INADINE® PVP-I Non-Adherent Dressing consists of a low adherent knitted viscose fabric impregnated with a polyethylene glycol (PEG) base containing 10% Povidone Iodine; equivalent to 1.0% available iodine.

Iodine is a highly effective topical antimicrobial that has been used clinically in the treatment of wounds for more than 170 years¹.

It has a broad spectrum of antimicrobial activity with efficacy against bacteria, microbacteria, fungi, protozoa and viruses¹.



PBI12176 